



Perfluoro-1-octanesulfonyl fluoride



Chemical Formula:	$\text{C}_8\text{F}_{18}\text{O}_2\text{S}$
Formula Weight:	502.12
Description:	Protecting Group & Tag
CAS Number:	307-35-7
Appearance:	Clear liquid
Soluble in:	Dichloromethane, chloroform, THF, ether, toluene and most other typical organic solvents
Stability:	Moisture sensitive

DESCRIPTION AND USES:


- Perfluorooctanesulfonyl fluoride is used to react with phenols to form perfluorooctanesulfonates.
- The perfluoroalkanesulfonyl tag has three functions:
 - A phenol protecting group
 - A phase tag for reaction mixture purification
 - A hydroxy activating group for the cross-coupling reaction
- Compounds bearing the tag can be manipulated in multi step synthesis. Detagging can be achieved by cross-coupling reactions to introduce additional diversity point or by traceless manner.²⁻⁵
- The intermediates can also be purified by fluoruous SPE¹ or other methods such as flash column chromatography.


TYPICAL TAGGING PROCEDURE:²

To a solution of a hydroxybenzaldehyde (11 mmol) in DMF (5 mL) was added K₂CO₃ powder (14 mmol) at room temperature. The reaction mixture was stirred for 10 minutes before the addition of perfluorooctanesulfonyl fluoride (9 mmol). The mixture was heated at 70°C for 2 hours. The purification of crude product can be conducted by F-SPE¹, precipitation, or flash chromatography on normal silica gel.

TYPICAL PROCEDURE FOR SUZUKI COUPLING:^{2,3}



A mixture of **2** (0.2 mmol), boronic acid (0.19 mmol), K₂CO₃ (0.4 mmol), Pd(dppf)Cl₂ (10 mol%) in 4:4:1 acetone/toluene/H₂O (1 mL) was heated at 130°C for 20 min under microwave irradiation (see figure below). The mixture was loaded onto a 5g FluoroFlash[®] SPE cartridge and then eluted with 80:20 MeOH/H₂O (10 mL). The MeOH/H₂O fraction was concentrated to afford the desired product **3**. The fluoruous component on the cartridge can be washed out by acetone (2x10 mL).



- i) 2-aminopyridine (1.1 equiv), cyclohexylisocyanide (1.2 equiv), Sc(OTf)₃ (0.05 equiv), 3:1 DCM/MeOH, microwave (150 w, 150 °C, 10 min), 76%;
- ii) 4-methoxyboronic acid (0.95 equiv), Pd(pddf)Cl₂ (0.1 equiv), K₂CO₃ (2.0 equiv), 4:4:1 acetone/toluene/H₂O, microwave (150 w, 130 °C, 20 min), 75%



TYPICAL PROCEDURE FOR THIOL COUPLING:^{2,4}

A mixture of **6** (0.2 mmol), thiol (0.19 mmol), K₂CO₃ (0.4 mmol), Pd(dppf)Cl₂ (10 mol%) in 4:4:1 acetone/toluene/H₂O (1 mL) was heated at 150°C for 10 min under microwave irradiation (see figure below). The mixture was loaded onto a 5g FluoroFlash[®] SPE cartridge and then eluted with 80:20 MeOH/H₂O (15 mL). The MeOH/H₂O fraction was concentrated to afford the desired product **7**.



TYPICAL PROCEDURE FOR TRACELESS DETAGGING:⁵

A mixture of **10** (0.2 mmol), HCOOH (5.8 mmol), K₂CO₃ (0.66 mmol), Pd(dppf)Cl₂ (13 mol%) in 4:4:1 acetone/toluene/H₂O (1 mL) was heated at 100°C for 20 min under microwave irradiation (see figure below). The mixture was loaded onto a 5g FluoroFlash[®] SPE cartridge and then eluted with 80:20 MeOH/H₂O (15 mL). The MeOH/H₂O fraction was concentrated to afford the desired product **11**.



REFERENCES:

- 1) FTI Application Note "Fluorous Solid Phase Extraction"
- 2) Lu, Y., Zhang, W., *QSAR & Combinatorial Chemistry*, in press.
- 3) Zhang, W., Chen, C. H.-T., Lu, Y., Nagashima, T., *Org. Lett.*, 2004, 6 (9), 1473.
- 4) Zhang, W., Lu, Y., Chen, C. H.-T., *Molecular Diversity*, 2004, 7, 199.
- 5) Zhang, W., Nagashima, T., Lu, Y., Chen, C., H.-T., *Tetrahedron Lett.*, 2004, 45, 4611.

FluoroFlash® is a registered trademark of Fluorous Technologies Inc.
20041026